


I'm not robot  reCAPTCHA

I'm not robot!

Sonics vibra cell manual

Sonics vibra cell vcx 130 manual. Vivo mobile vibration not working. Sonics vibra-cell vcx 750 manual. Sonic vibration definition. Sonics vibra cell vc 505 manual.

LABORATORY SCALE MODELS (VOLUME 0.2mL-1L*) LABORATORY SCALE MODELS (VOLUME 0.2mL-1L*) LABORATORY SCALE MODELS (VOLUME 0.2mL-1L*) The VCX 750 is a powerful and versatile ultrasonic liquid processor. Using a variety of accessories, the system can process a wide range of sample types and volumes for many different applications. The VCX 750 is microprocessor based and programmable.



Intuitive screens are user friendly and easy to navigate.



The 750 watt power supply has the capability of operating our largest assortment of horns and accessories. A standard 1/2" diameter probe (50 - 250mL volumes) with replaceable tip is included with the system. Part No. VCX 750A can be used to order the system without the standard probe. The VCX 500 offers 500 watts of power to handle a range of liquid processing applications.



Artisan Technology Group

This model shares the same microprocessor based operating system and programming screens as the VCX750 system. Using a variety of accessories, the system can process a wide range of sample types and volumes for many different applications. The VCX 500 includes a standard 1/2" diameter probe (50 - 250mL volumes) with replaceable tip. Part No. VCX 500A can be used to order the system without the standard probe. The VC 505 is a cost effective alternative to the VCX model. This 500 watt system offers the same processing power and choice of accessories as the VCX500 without temperature monitoring or energy set point programming. A standard 1/2" diameter probe (50 - 250mL volumes) with replaceable tip is included with the VC 505. Many optional accessories are available to meet the needs of any application and additional sample volumes. Part No. VC 505A can be used to order the system without the standard probe.



UNIQUE SONIC VIBRATION EXERCISE & THERAPY Sonix Whole Body Vibration machines use advanced sonic science to deliver precise waves of infrasound and super-low frequencies to the entire body via electromagnetic technology. The principle of sonic vibration makes Sonix units the world's first devices optimized for vertical vibration motion, and resolves the problems with noise and component wear of other motorized technologies. Designed for professional or home use, the PROFESSIONAL model has 12 Automatic Modes and a Manual Mode, with a weight capacity of up to 551 lbs. Its Dual Mode Operation audio system allows the user to feel the selected vibration frequency as well as the sound from external sources at the same time. Audio jacks permit sound file and music input to the Sonix unit from MP3/iPod players, PCs and frequency generators (e.g. for RIFE frequencies). WHAT IS PEMF (Pulsed Electromagnetic Field)? PEMF devices are used worldwide to accelerate healing and reduce pain or discomfort. Only Sonix machines combine the power of PEMF with sonic vibration therapy. PEMF alkalizes the blood and improves blood flow while increasing oxygen uptake. Most good PEMF units start at \$6000 and are about 1/40th the power of the Sonix. See About PEMF for further information. THE FANTASTIC FEELING OF SONIC WAVES Sonic vibration technology provides enjoyable sensations from the fantastic feeling of Sonic Waves passing through the body. The therapeutic sound frequencies are highly effective in providing both a full-body workout and intensive exercise for selected body parts while in an upright position. Frequencies (cycles per second) are user selectable from 3-70Hz with adjustable intensity. MAXIMUM EFFECT IN MINIMUM TIME Sonic wave exercise with user-selected frequency, intensity and duration can increase maximum strength by over 30% when combined with an existing training program. Traditional exercise time can be reduced by up to 85%, and sonic vibration therapy can help to prevent injuries that may occur during weight training. A 10-minute Sonix exercise program can provide a highly effective workout. SONIC EXERCISE BENEFITS THE ENTIRE BODY Defined sonic wave exercises stimulate fine muscle as well as large muscles, including muscle groups and body areas that are not easily exercised in traditional programs for fast stretching and recovery from muscle fatigue. Sonic Wave therapy can stimulate the facial muscles, bone and tissue cells, as well as the whole body to provide the most unique and thorough workout you've ever experienced. EASY, SAFE & EFFECTIVE FOR EVERYONE Sonic Wave vibration technology makes highly effective exercises available without burdening joints or ligaments. In a small space, safe programs can be easily developed and utilized from children to seniors, by the healthy or those with limited mobility, and for the benefit of deconditioned and super athletes alike. CUSTOMIZED PROGRAMS FOR INDIVIDUAL NEEDS Depending on the user's physical condition and any medical or other considerations, the intensity and duration of a selected sonic vibration frequency can be finely controlled and customized for a self-paced movement program. It can be applied to a wide variety of desired outcomes in diverse fields including health and wellness, beauty, medical issues, general fitness and professional sports. SONIX PROFESSIONAL VC 15 SPECIFICATIONS Size 870 x 910 x 1500 (mm) / 34.25" x 35.8" x 59" Weight: 85Kg / 187.4 lbs. Weight limit: 250Kg / 551 lbs. Frequency Range: 3-70Hz Power source: 110-120V 50/60Hz Power consumption: 400W Operating environment: Temperature 0-40°C (32-104°F); Humidity : 0-80% Sonic wave generator: 1500 x 2 Programs: 12 pre-programmed modes; 1 manual mode Intensity: Adjustable Assembly time: approx. 7 minutes (2 people) Color: Silver Warranty: 3 years labor. 1 year parts. NOTE: free delivery is "curbside delivery." Two boxes on a pallet to your driveway. Driver not required to assist with moving boxes. The manual will be emailed for usage. © 2023 Sonics & Materials, Inc.